

To Whom it May Concern:

All Biomeme Inc products have been cleared for international travel.

The Harmonized Tariff Schedule codes for our products are as follows:

- Franklin Thermocycler: 9027.80.90
- Android Smartphone: 8517.12.00
- Lab Reagents: 3822.00.50.90

Furthermore, our products meet the following harmonized standards:

Safety:	IEC 61010-1:2010, AMD1:2016 IEC 61010-2-010:2014 IEC 61010-2-081:2015 IEC 62133-2:2017
EMC:	IEC 62479:2010 IEC 61326-1:2013 ETSI EN 301 489-1 V2.1.1 ETSI EN 301 489-17 V3.1.1
Hazardous Substances:	EN 50581:2012

Battery Type	UN ID	Watt-Hour Rating	Battery Weight (kg)	Battery Location	Shipment Type
Lithium Ion Battery	UN 3481	96.5 Wh	0.410	Contained in Equipment	Excepted

Therefore:  $5.0/0.410 =$  **packages of up to 12 Franklin devices may be transported via passenger air/rail.**

Please contact us at [support@biomeme.com](mailto:support@biomeme.com) or call +1 267-930-7707 if further information is needed.  
Thank you!